



# Cleaning, Disinfecting and Sterilization Resistance of 3M™ Medical Assembly Vinyl Tape 471 Solution for Reusable Medical Devices

Technical Bulletin

January, 2013

## Introduction:

With the development of reusable medical devices, sterilization between uses has become an important consideration in the development of medical devices. In May 2011, the FDA published a Draft Guidance for Industry and FDA staff entitled *Processing/Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling*. When finalized, medical device manufacturers will need to include validated labeling instructions for reprocessing reusable medical devices in premarket submissions.

3M™ Vinyl Tape 471 is used in the manufacture of medical devices. This tape is “colored throughout” which resists scrapes, chipping and cracking through repeated sterilization cycles as well as cleaning and disinfecting solutions. It is easily applied color coding, to help organize and track surgical instruments. It is a transparent pressure sensitive synthetic rubber adhesive for high initial adhesion. With the shift to reusable medical devices, these devices must be able to withstand the specified sterilization methods outlined by industry requirements. 3M has tested the integrity of our product by conducting evaluations to several industry protocols.

## Protocol Description:

To meet the requirements of various types of medical devices different cleaning, disinfecting, and sterilization techniques were tested including heat, radiation and chemical immersion. The tests conducted were 90° Angle Peel Adhesion tests based on ASTM D3330 on stainless steel panels. Upon completion of the test methods, the samples were visually inspected after each cycle for edge lift and buckling. All tests were conducted on the exposed samples after a minimum of 24 hours at ambient temperature. The different protocols used are described below:

| Sterilization Method            | Description   | Protocol  | Visual Inspection |
|---------------------------------|---|---|-------------------|
| Chemical Disinfecting Solutions | <ul style="list-style-type: none"><li>• 3.4% Gluteraldehyde Solution</li><li>• 0.55% o-Phthaldehyde Solution</li><li>• 1% Hydrogen Peroxide and 0.8% Peracetic Acid Solution.</li></ul> | For each solution, samples were immersed in the solution for 1 hour for 5 cycles, followed by a rinse with distilled water. | All Pass          |
|                                 | <ul style="list-style-type: none"><li>• Quaternary Ammonium Chloride Solution</li></ul>   | Immersed for 1 minute for 5 cycles, followed by a rinse with distilled water.   | Pass              |
| Steam Autoclave                 | <ul style="list-style-type: none"><li>• Steam under pressure</li><li>• Temperatures between 100-135°C</li></ul>   | Exposed to 3 cycles @ 250°F (121°C) for 20 minutes under 25 psi (172 kPa) pressure.   | Pass              |
| Gamma                           | <ul style="list-style-type: none"><li>• Gamma rays</li><li>• Cobalt 60 radiation source</li><li>• Ambient temperature</li></ul>   | Sample was exposed to 1 cycle 25.0-45.0 kGy   | Pass              |
| Ethylene Oxide                  | <b>Gas Exposure</b> <ul style="list-style-type: none"><li>• Temperature range between 86°F (30°C) and 140°F (60°C)</li><li>• Gas concentration between 200 and 800 mg/l</li></ul>       | Sample was exposed to 1 cycle at 131°F (55°C) for 60 minutes  | Pass              |


**Note:** The following technical information and data should be considered representative or typical only and should not be used for specification purposes.

# Technical Bulletin

## Cleaning, Disinfecting and Sterilization Resistance of 3M™ Medical Assembly Vinyl Tape 471 Solution for Reusable Medical Devices

### Results:

For the vinyl tape product, the results for all cycles were averaged. The mean value was compared to the corresponding control value and the percent difference is presented below:



In these evaluations, performance of the tape was not negatively affected by exposure to the various sterilization methods. The product performed similar to the control sample, with all samples maintaining at least 85% of the strength of the control.

# Technical Bulletin

## Cleaning, Disinfecting and Sterilization Resistance of 3M™ Medical Assembly Vinyl Tape 471 Solution for Reusable Medical Devices

---

### Technical Information

The technical information, recommendations and other statements contained in this document are based upon tests or experience that 3M believes are reliable, but the accuracy or completeness of such information is not guaranteed.

---

### Product Use

Many factors beyond 3M's control and uniquely within user's knowledge and control can affect the use and performance of a 3M product in a particular application. Given the variety of factors that can affect the use and performance of a 3M product, user is solely responsible for evaluating the 3M product and determining whether it is fit for a particular purpose and suitable for user's method of application.

---

### Warranty, Limited Remedy, and Disclaimer

Unless an additional warranty is specifically stated on the applicable 3M product packaging or product literature, 3M warrants that each 3M product meets the applicable 3M product specification at the time 3M ships the product. 3M MAKES NO OTHER WARRANTIES OR CONDITIONS, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OR CONDITION OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY IMPLIED WARRANTY OR CONDITION ARISING OUT OF A COURSE OF DEALING, CUSTOM OR USAGE OF TRADE. If the 3M product does not conform to this warranty, then the sole and exclusive remedy is, at 3M's option, replacement of the 3M product or refund of the purchase price.

---

### Limitation of Liability

Except where prohibited by law, 3M will not be liable for any loss or damage arising from the 3M product, whether direct, indirect, special, incidental or consequential, regardless of the legal theory asserted, including warranty, contract, negligence or strict liability.



This Industrial Adhesives and Tapes Division product was manufactured under a 3M quality system registered to ISO 9001 standards.



### Industrial Adhesives and Tapes Division Converter Markets

1030 Lake Road  
Medina, OH 44256-0428  
800-422-8116 • 877-722-5072 (fax)  
www.3M.com/converter



*Recycled Paper*  
40% pre-consumer  
10% post-consumer

3M is a trademark of 3M Company.  
Printed in U.S.A.  
©3M 2013 (1/13)